E sanabelle HME-A6

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GB sanabelle HME-A6

Artificial nose with safety valve and oxygen connector

Indication

For the efficient filtering, warming and humidifying of inhaled air for patients with spontaneous breathing through an endotracheal tube or tracheal cannula, in hospital or at home.

Description

The sanabelle HME-A6 comprises a plastic housing with an insert of hydroscopic foam material. The HME has a 15 mm female connector, which can be fitted to the standardised 15 mm connector of the tracheal cannula. The HME is latex-free, DEHP-free and electrically non-conductive.

The safety valve minimizes the risk of occlusion of the humidifying medium, allows short-term high tidal volume and coughing up of larger quantities of secretions.

Usage

Unpack the sanabelle HME-A6 only before immediate use carefully. Fit the HME to a tracheal cannula with a 15 mm connector. The casing is to be fitted to the open end of the endotracheal tube or tracheostomy cannula. Make sure the connection is securely attached! If supplemental oxygen has to be delivered, connect the oxygen connector to an oxygen supply. The normal humidifying effect sets in after only a few breaths.

The HME must not be used together with a humidifier operating with heated or nebulised water or in combination with a nebuliser.

IMPORTANT: The instructions for use of the sanabelle HME-A6 must be carefully read prior to use, to minimise the potential risk of dyspnoea, labored breathing, hypercapnia and hypoxia.

Contraindications

- The sanabelle HME-A6 must not be used outside the recommended tidal volume range, in order to prevent an increased concentration of CO₂ in case of a low tidal volume, or insufficient humidification, if the tidal volume is too high. Please observe the doctor's recommendation. Please consult a specialist in case of questions.
- The device is not to be used on patients with strong exsiccosis or excessive secretions in the lungs and airways or in case of leaking gas flow (expired gas does not flow through the humidifier).
- Do not connect the HME to an artificial ventilation system.

Precautions

- Always ensure the patency of the HME.
- Replace the sanabelle HME-A6 after a maximum of 24 hours to counteract any possible bacterial or fungicidal growth. The HME must be replaced more frequently depending on the needs, for example, if the breathing resistance is in danger of increasing due to accumulations of secretion.
- Monitor patients for potential fatigue! Upon indication of difficulty breathing, remove the HME immediately! Reason: the HME's resistance requires increased breathing effort.
- The sanabelle HME-A6 is for single use only! Cleaning, disinfecting or re-sterilising it for reuse may impair the function and increase the risk of transmitting infections.
- Dispose of used HMEs immediately in a hygienically proper manner!

Note

It is recommended to assure the supervision of children and people with intellectual disabilities in order to prevent the risk of an accidental decannulation.

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