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
GB sanabelle TM100 cuff pressure monitor

Intended Use

The sanabelle TM100 cuff pressure monitor is used for ventilation and pressure control of tracheal cannulae and endotracheal tubes with low pressure cuffs by health professionals or users having received respective training.

Warning

The cuff pressure monitor is intended for use with tracheal cannulae or endotracheal tubes. The instructions for use and warnings of the product used with the cuff pressure monitor must be considered, too.

The cuff pressure monitor is not suitable for MRI.

For cleaning adhere to the respective cleaning instructions!


In case that the cuff pressure monitor is used for monitoring, i.e. remains connected to the tracheal cannula or endotracheal tube, the pressure must be checked regularly every 15 minutes. If necessary pressure can be adjusted by use of the knurled screw or the inflation bulb as the cuff pressure monitor cannot hold the cuff pressure constantly.

To prevent damage it has to be observed that the Luer connector of the cuff pressure monitor is not subject to prying by the connection to the cannula, endotracheal tube or the connecting tube.

If the cuff pressure monitor is used with a connecting tube without its own one-way-valve, it has to be disconnected together with the connecting tube as otherwise the system is open.

If a connecting tube is used, it has to be observed that it is not buckled during its use. The sanabelle TM100 cuff pressure monitor is intended for use by health professionals or users having received respective training.



Testing before each use

Before each use the cuff pressure monitor has to be tested as follows:

1. The needle must not deflect more than two graduation marks from the zero position. If a greater deviation is determined, the (decontaminated) device has to be sent in to medi1one for repair.
2. Close the release valve by turning the knurled screw clockwise.
3. Seal the Luer connector with a finger minding proper hygiene. (Figure 2)
4. Inflate with sealed Luer connector to a pressure of approx. 40 cmH₂O with the inflation bulb. This pressure must be held constant for 2-3 seconds. If the pressure decreases, the (decontaminated) device has to be sent in to medi1one for repair.

If a connecting tube is used, the following additional tests have to be performed:

5. Attach the connecting tube to the Luer connector of the cuff pressure monitor. (Figure 3)

6. Seal the end of the connecting tube with a finger. (Figure 4)

7. Inflate with sealed Luer connector to a pressure of approx. 40 cmH₂O with the inflation bulb. This pressure must be held constant for 2-3 seconds. If the pressure decreases, the connecting tube may be leaky.

In the latter case the test has to be performed as previously described without the connecting tube to rule out the cuff pressure monitor as error source. If this test confirms that the connecting tube is leaky, the connecting tube has to be replaced and the test has to be repeated with the new connecting tube.

Instructions for use

Before use mind the instructions for use of the respective tracheal cannula or endotracheal tube.

1. For use close the release valve by turning the knurled screw (Figure 5) and connect the cuff pressure monitor – where applicable with a connecting tube – to the pilot balloon of the tracheal cannula or the endotracheal tube. It has to be observed for all connections with a one-way-valve that the valve is opened, because with closed valve the cuff pressure will not be displayed correctly.

2. If not indicated differently in the instructions for use of the respective tracheal cannula or endotracheal tube, the cuff should be inflated to a pressure of 60 to 90 cmH₂O after intubation/insertion to ensure that it is in close contact with the tracheal wall.
3. Release the pressure immediately by turning the knurled screw (Figure 5) until the cuff pressure is between 20 cmH₂O and 30 cmH₂O (dark green/graphically offset area with labelling on the scale). (Figure 1)

This pressure range is chosen to reduce on the one hand the risk of aspiration or pneumonia by too low pressure and on the other hand the risk of ischemia of the tracheal mucous membrane (e.g. Hamilton & Grap, 2012 Heart Lung 41: 167-72).

Attention: The values provided are guidance values and have to be adjusted to the recommendations of the attending physician and/or the individual circumstances. E.g. it can be necessary for the use with children to adjust to lower pressure values.

Because the cuff pressure may change over time or after repositioning of the patient, periodic controls of pressure approximately every four hours are recommended. The cuff pressure can be reduced by opening the release valve with the knurled screw or increased by operating the inflation bulb if necessary during controls.

Before extubation the cuff of the tracheal cannula or endotracheal tube should be deflated completely with a syringe according to the respective instruction for use by the manufacturer.

The measurement accuracy of the cuff pressure monitor is ± 2 cmH₂O for the whole measuring range.

Remark on scaling:

The following formula applies to the conversion of cmH₂O into the SI unit Pascal:

$$1 \text{ cmH}_2\text{O} = 0.98067 \text{ hPa}$$

Cleaning instructions

If necessary the device can be cleaned with a mild soap solution.

No aggressive cleaning agents must be used. The device must not be immersed in liquids. For disinfection the device can be wiped with a commercially available disinfectant on alcohol base. The cuff pressure monitor cannot be automatically re-processed or sterilized.

Maintenance

Maintenance of the device by the user is not intended.

If deviations are detected during the tests or during use that exceed the given

tolerance values, the decontaminated product must be sent in to medi1one for repair.

Storage

The product must be stored cool and dry. Protect from sunlight and extreme cold. Short term storage for a maximum of 30 days at a maximum temperature of 50°C is possible.

Disposal

The product can be disposed after decontamination in accordance with the valid legal regulations.

Ordering information

REF	Article	Sales unit
77701100	sanabelle TM100 cuff pressure monitor	1 piece
77701199	sanabelle connecting tube for TM100	1 piece or 10 pieces
77701197	Protective ring (Spare part)	1 piece

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